

Agilent's Solutions for Nanoparticles Analysis

ICP-MS, ICP-MS/MS, software



Alain Desprez
Agilent Technologies



The Total Solution

Choice of Mass Spectrometers

- 7800 quadrupole ICP-MS – cost-effective, 3 ms dwell time
- 7900 quadrupole ICP-MS – ultra-sensitive, ultra-fast 0.1 ms dwell time
- 8900 ICP-QQQ – highest sensitivity and lowest background, MS/MS for interference-free determination of SiO₂ and TiO₂ nanoparticles

Fully automated acquisition and data analysis configuration within MassHunter 4.3 software guided by Method Wizard

- Supports single particle mode with 1 or 2 elements and field-flow fractionation (FFF) coupling



7800 ICP-MS




7900 ICP-MS



8900 ICP-QQQ

Single Particle-ICP-MS (spICP-MS)

- Uses normal sample intro system and TRA data collection method
- If sample is sufficiently dilute, each nanoparticle gives a distinct transient signal
- Particle concentration, particle size (diameter), size distribution and composition are obtained



Method Wizard now supports nano analysis by both sNP and FFF modes

MassHunter Method Wizard

Analysis Mode

Select appropriate Analysis Mode. Also make the following changes to radio button selections.

Conventional Analysis

- Spectrum/TRA
Spectrum: Conventional Mass Spectra are acquired.
TRA: Conventional Time Resolved Data are acquired.

Nanoparticle

- Single Particle Analysis
Single Particle Data are acquired using fast TRA mode.
- Multiple Elements within Single Particle
2 elements within a single particle are acquired using fast TRA mode.
- FFF
Particles are acquired using Field-Flow Fractionation.

To continue, click Next.

[Help](#)

Method Wizard now supports nano analysis by both sNP and FFF modes

Most parameters are automatically entered after providing a couple of simple inputs such as sample tube ID and Reference Material Name

MassHunter Method Wizard

Analysis Mode

Select appropriate Analysis Mode. You can make the following changes to radio button selections.

Conventional Analysis

- Spectrum/TRA
Spectrum: Conventional
TRA: Conventional

Nanoparticle

- Single Particle Analysis
Single Particle Data mode.
- Multiple Elements
2 elements within a single particle are acquired using fast scan.
- FFF
Particles are acquired using Flow Fractionation.

Single Particle Analysis Configuration

Set parameters for Single Particle Analysis.

Sample Pump Tube ID: 1.02 mm

Sample Inlet Flow: 0.346 ml/min

Response Factor Calibration Solution:

Response at 107 amu: 200000 cps/ppb

Reference Material: NIST RM 8012

Reference Element Mass: 197 amu

Mean Reference Particle Diameter: 30 nm

Reference Element Density: 19.32 g/cm³

Mass Concentration of Reference Material: 5.0 ng/l


Unknown Sample:

Target Element Mass: 107 amu

Analyte Mass Fraction: 1.000

Analyte Element Density: 10.50 g/cm³

Help < Back Next > Finish Cancel



Just Complete the Sample List, Add it to the Queue...


Skip	Sample Type	Sample Name	Comment	Vial#	Replicates	Dilution
<input type="checkbox"/>	IonicBlk	Ionic Blank				
<input type="checkbox"/>	IonicStd	Ionic Standard	10 ppt Ag			
<input type="checkbox"/>	RM	NIST RM 8012	30nm RM			
<input type="checkbox"/>	Sample	Unknown 1				
<input type="checkbox"/>	Sample	Unknown 2				
<input type="checkbox"/>	Sample	Unknown 3				
<input type="checkbox"/>	Sample	Unknown 4				
<input type="checkbox"/>	Sample	Unknown 5				
<input type="checkbox"/>	Sample	NIST RM 8012				
<input type="checkbox"/>						

And Review the Results

Sample	Acq. Date-Time	Type	Sample Name	Nebulization Efficiency	# of Particles	Conc. (particles)	Conc. (ng/l)	Ionic Conc. (ppb)	BED (nm)	Particle Size (nm)
1	10/29/2014 12:39:00 AM	RM	Au 60nm 100ppt	0.059	1173	5.6E+7	100.0	0.0645	45.12	55
2	10/29/2014 12:30:19 AM	Sample	Au 60nm 100ppt	0.059	758	4.5E+7	99.9	0.0683	39.43	59
3	10/29/2014 12:34:18 AM	Sample	Au 60nm 100ppt	0.059	954	4.6E+7	90.0	0.0519	31.91	57
4	10/29/2014 12:41:44 AM	Sample	Au 60nm 100ppt	0.059	1207	5.8E+7	110.5	0.0611	39.26	56
5	10/29/2014 12:44:24 AM	Sample	Au 60nm 100ppt	0.059	1207	5.8E+7	110.4	0.0596	34.12	56
6	10/29/2014 12:47:16 AM	Sample	Au 60nm 100ppt	0.059	1188	5.7E+7	108.9	0.0593	31.05	56
7	10/29/2014 12:49:54 AM	Sample	Au 60nm 100ppt	0.059	1223	5.9E+7	109.7	0.0598	35.77	56
8	10/29/2014 12:52:23 AM	Sample	Au 60nm 100ppt	0.059	1141	5.5E+7	104.2	0.0581	30.76	56

100 μ s integration time and no settling time with the Agilent 8900 and 7900 captures all the data


Zoom in...



Raw Data to Right Results in Seconds



Composite Sample	Observed concentration (particles/L)	Observed concentration (ng/L)	Prepared total concentration (ng/L)	Recovery
60nm (NIST 8013) 50 ng/L + 30nm (NIST 8012) 5 ng/L	4.78×10^7	57.6	55	105%



Reference particle size by TEM (nm)	Observed particle size (nm)	Prepared particle number ratio (%)	Observed particle number ratio (%)
56.0±0.5	58	55.6	57.8
27.6±2.1	28	44.4	42.2

Determining Really Small Nanoparticles

What are the limitations and how can ICP-MS/MS solve them?

There are only 2 things that limit the minimum size nanoparticle that can be detected in single particle mode.

1. Sensitivity
2. Background (raw instrument background, spectral interferences, dissolved ionic background)

Everything else is secondary.

Nanoparticle size (diameter) is a function of the cube root of the mass. If the diameter of a particle is reduced by half, the mass and therefore the signal is reduced by 8X. HIGH SENSITIVITY IS THE CRITICAL FACTOR FOR DETECTING SMALL PARTICLES.

Background is also critical because nanoparticle peaks are detected as x times the signal to noise. Small peaks can be lost in the noise.



Agilent application note 5991-6596

Challenge of Si or Ti Measurement by Conventional (Single) Quadrupole ICP-MS

- **Interferences**

- Si is interfered by C, N and O overlaps (ubiquitous in most aqueous/organic samples)
- Ti is interfered by polyatomic overlaps in P, S, Si and C matrices

These interferences make the analysis of SiO₂ or TiO₂ NPs in real samples difficult


- **⁴⁸Ti (74% abundance) cannot be measured reliably by ICP-QMS, due to the ⁴⁸Ca isobaric interference**

- Less interfered ⁴⁹Ti (5.4% abundance) has less sensitivity, so detection of smaller particles is compromised

➔ ICP-MS/MS can solve these problems


ICP-MS/MS Technique for Si (On-Mass Method)

ICP-MS/MS uses two quadrupoles separated by a collision/reaction cell




SiO₂ Nanoparticle Reference Material Results

50 nm




Nominal size (nm)	Median size (nm)	Most Freq. Size (nm)	Mean Size (nm)	TEM Diameter (1) (nm)	BED (2) (nm)
50	49	50	49	46.3 ± 3.1	23
60	59	58	60	57.8 ± 3.5	22
100	99	100	100	97.0 ± 4.8	25
200	199	204	199	198.5 ± 10.5	25


60 nm



100 nm





200 nm



- (1) values supplied by nanoComposix
- (2) Background equivalent diameter ~ size detection limit – most likely limited by dissolved Si.

ICP-MS/MS Technique for Ti (Mass-Shift Method)



TiO₂ analysis - Interference Removal by O₂-H₂ Mode

Determining interferences in ionic matrix solutions

Mode	Sensitivity (cps/ppb)	BEC (ppb)	Quantitative result as Ti (ppb)					
			P 100 ppm	S 100 ppm	Ca 50 ppm	Si 50 ppm	Ethanol 0.1%	Matrix mixture
No gas	155000	0.016	1.7	6.0	225	0.39	0.14	261
O ₂ + H ₂	79000	0.001	0.010	0.001	0.18	0.054	0.001	0.023

No gas : Single Quad mode, Q2=48

O₂ + H₂ : MS/MS mode, Q1=48, Q2=64


Advantage of O₂ + H₂ mode

- Can use the most abundant isotope of Ti (mass 48, 73.7% relative abundance) for highest sensitivity
- Can remove almost all interferences including ⁴⁸Ca

TiO₂ NPs in Matrix Samples


O₂-H₂ MS/MS mode

Sunscreen
in DI water




- Mean size : 77 nm
- **BED : 16 nm**

Sunscreen
in tap water



- Mean size : 79 nm
- **BED : 10 nm**

Sunscreen
in matrix mixture
(100 ppm S & P,
50 ppm Ca & Si,
0.1% ethanol)



- Mean size : 84 nm
- **BED : 22 nm**

Conclusions

The sp-ICP-MS/MS method with MS/MS for SiO₂ and TiO₂ nanoparticle analysis:

- Provides high sensitivity for Si and Ti
- Effectively eliminates polyatomic ions which interfere with Si or Ti analysis thereby significantly lowering the minimum detectable particle size even in complex matrices
- Permits the use of major isotopes (²⁸Si, ⁴⁸Ti), even in the presence of polyatomic and isobaric (⁴⁸Ca) overlaps
- Provides quick and accurate results for SiO₂ and TiO₂ particles smaller than 100 nm

Thank You